

SELE / CD62E / E-selectin Mouse anti-Human Monoclonal (FITC) Antibody - LS-C44155 - LSBio	
<b>CatalogID:</b>	LS-C44155
<b>Target:</b>	selectin E
<b>Synonyms:</b>	SELE Antibody, CD62E Antibody, ELAM Antibody, ESEL Antibody, E-selectin Antibody, LECAM2 Antibody, CD62E antigen Antibody, ELAM-1 Antibody, ELAM1 Antibody, Selectin E Antibody
<b>Host</b>	SELE antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	SELE / CD62E / E-selectin antibody was raised against Human
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Tris-sodium chloride, 0.02% sodium azide, 1% BSA.
<b>Usage Summary:</b>	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100 ul. Method sheets are available on request.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	