

IL2RA / CD25 Mouse anti-Human Monoclonal (RPE) Antibody - LS-C44320 - LSBio

CatalogID:	LS-C44320
Target:	interleukin 2 receptor, alpha (IL2RA)
Synonyms:	IL2RA Antibody, CD25 antigen Antibody, CD25 Antibody, IDDM10 Antibody, IL-2 receptor subunit alpha Antibody, IL2R Antibody, Interleukin 2 receptor, alpha Antibody, IL-2-RA Antibody, IL2-RA Antibody, TCGFR Antibody, IL-2R subunit alpha Antibody, p55 Antibody, TAC antigen Antibody
Family / Subfamily:	Interleukin
Host	IL2RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	IL2RA / CD25 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	IL2RA / CD25 antibody was raised against human PHA blasts, day 3 of culture.
Specificity:	Recognizes the human CD25 cell surface antigen, a 55kD glycoprotein which is a low affinity receptor for IL-2. CD25 is expressed by activated T lymphocytes
Reactivity:	Human
Purification:	Protein A purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, PBS, pH 7.4, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences