

CD4 Mouse anti-Chicken Monoclonal (RPE) Antibody - LS-C44343 - LSBio

CatalogID:	LS-C44343
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD4 antibody was raised against Chicken
Antigen Type:	Cells
Immunogen:	CD4 antibody was raised against chicken embryonic thymocytes.
Specificity:	Recognizes the chicken homologue of human CD4, a 64kD cell surface protein expressed by thymocytes and a subset of T cells
Reactivity:	Chicken, Turkey
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1 - 1:5) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences