

CD4 Mouse anti-Sheep Monoclonal (RPE) Antibody - LS-C44424 - LSBio

CatalogID:	LS-C44424
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD4 antibody was raised against Sheep
Immunogen:	CD4 antibody was raised against fetal thymocytes.
Specificity:	Recognizes the ovine CD4 cell surface glycoprotein, which is expressed by a subset of mature T lymphocytes
Reactivity:	Sheep, Goat
Purification:	Protein G purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, PBS, pH 7.4, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1 - 1:5) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences