

CD45 Mouse anti-Human Monoclonal (APC) (MEM-28) Antibody - LS-C46194 - LSBio

CatalogID:	LS-C46194
Target:	protein tyrosine phosphatase, receptor type, C (PTPRC)
Synonyms:	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
Family / Subfamily:	Protein Phosphatase / R1/R6
Host	PTPRC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-28
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	CD45 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45 antibody was raised against human thymocytes and T lymphocytes.
Specificity:	Reacts with all alternative forms of human CD45 antigen (Leukocyte Common Antigen), a 180-220 kD single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except erythrocytes and platelets.
Reactivity:	Human
Purification:	Ascites
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences