

**CD45RO Mouse anti-Human Monoclonal (APC) (UCHL1) Antibody - LS-C46254 - LSBio**

<b>CatalogID:</b>	LS-C46254
<b>Target:</b>	CD45RO
<b>Host</b>	CD45RO antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	UCHL1
<b>Conjugations:</b>	Allophycocyanin (APC)
<b>Immunogen Species:</b>	CD45RO antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45RO antibody was raised against human IL-2 dependent T cells.
<b>Specificity:</b>	Recognizes CD45RO, a 180 kD low molecular weight isoform of the leukocyte common antigen (LCA). The antigen is expressed on a subset of memory/activated T cells and on cortical thymocytes.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences