

CD7 Mouse anti-Human Monoclonal (FITC) (MEM-186) Antibody - LS-C46340 - LSBio

CatalogID:	LS-C46340
Target:	CD7 molecule
Synonyms:	CD7 Antibody, CD7 antigen Antibody, CD7 antigen (p41) Antibody, CD7 molecule Antibody, LEU-9 Antibody, T-cell leukemia antigen Antibody, T-cell antigen CD7 Antibody, T-cell surface antigen Leu-9 Antibody, TP41 Antibody, gp40 Antibody, p41 protein Antibody, Tp40 Antibody
Host	CD7 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	MEM-186
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD7 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD7 antibody was raised against human acute myelogenous leukaemia cell line KG-1.
Specificity:	The MEM-186 antibody reacts with CD7, a 40 kD type I transmembrane glycoprotein expressed on peripheral blood T lymphocytes, NK-cells, hematopoietic progenitors, monocytes (weakly) and also on acute lymphocytic leukemia
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences