

GYPA / CD235a / Glycophorin A Mouse anti-Human Monoclonal (N-Terminus) (FITC) (HIR2) Antibody - LS-C46406 - LSBio	
CatalogID:	LS-C46406
Target:	glycophorin A (MNS blood group) (GYPA)
Synonyms:	GYPA Antibody, CD235a Antibody, CD235a antigen Antibody, Glycophorin A (MN blood group) Antibody, Glycophorin A, GPA Antibody, Glycophorin-A Antibody, GPA Antibody, GPERik Antibody, GPSAT Antibody, Glycophorin MiV Antibody, Glycophorin SAT Antibody, HGpMiX Antibody, HGpMiXI Antibody, HGpMiV Antibody, HGpSta(C) Antibody, Glycophorin Mil Antibody, GpMiIII Antibody, MN Antibody, HGpMiIII Antibody, MN sialoglycoprotein Antibody, Sialoglycoprotein alpha Antibody, Mi.V glycoprotein (24 AA) Antibody, PAS-2 Antibody, Glycophorin A Antibody, Glycophorin Sta type C Antibody, MNS Antibody
Host	GYPA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	HIR2
Conjugations:	Fluorescein (FITC)
Immunogen Species:	GYPA / CD235a / Glycophorin A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GYPA / CD235a / Glycophorin A antibody was raised against synthetic peptide (Human, N-terminal).
Specificity:	Recognizes N-terminal, homologous portion of glycophorins A (GPA) and B (GPB) which are single-pass membrane sialoglycoproteins expressed on red blood cells and erythroid precursor cells. GPA and GPB provide the cells with a large mucin-like surface. The antibody HIR2 agglutinates untreated RBCs but it fails to agglutinate papain-treated cells. The antibody HIR2 significantly binds to GPA, but weakly to GPB. The antibody is useful in erythroid cell development studies, because HIR2 antigen is expressed on early erythroblasts, late erythroblasts, erythroblasts, mature erythrocytes and the cells of erythroid cell lines K562 and HEL, but not on all other cells (mature erythrocytes are characteristically CD235a positive and CD45 and CD71 negative)
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Ascites
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences