

CD45RO Mouse anti-Human Monoclonal (FITC) (UCHL1) Antibody - LS-C46422 - LSBio

CatalogID:	LS-C46422
Target:	CD45RO
Host	CD45RO antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	UCHL1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45RO antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RO antibody was raised against human IL-2 dependent T cells.
Specificity:	Recognizes CD45RO, a 180 kD low molecular weight isoform of the leukocyte common antigen (LCA). The antigen is expressed on a subset of memory/activated T cells and on cortical thymocytes.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences