

Human Kappa Light Chain Mouse anti-Human Monoclonal (FITC) (A8B5) Antibody - LS-C46424 - LSBio	
CatalogID:	LS-C46424
Target:	Human Ig Kappa Light Chain
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	A8B5
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Human
Specificity:	Reacts with kappa light chains (22.5 kD) of immunoglobulins.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	