

CD3 Mouse anti-Human Monoclonal (PE) (MEM-57) Antibody - LS-C46430 - LSBio

CatalogID:	LS-C46430
Target:	CD3
Host	CD3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MEM-57
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD3 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD3 antibody was raised against human thymocytes and T lymphocytes.
Specificity:	Reacts with gamma-epsilon and delta-epsilon dimers of human CD3 complex, a part of a bigger multisubunit T cell receptor complex (CD3/TCR) expressed on peripheral blood T lymphocytes and mature thymocytes
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences