

LFA-1 / CD11A Mouse anti-Human Monoclonal (PE) (MEM-25) Antibody - LS-C46437 - LSBio	
CatalogID:	LS-C46437
Target:	integrin, alpha L (antigen CD11A (p180), lymphocyte function-associated antigen 1; alpha polypeptide) (ITGAL)
Synonyms:	ITGAL Antibody, CD11A Antibody, Integrin alpha-L Antibody, LFA-1 alpha Antibody, LFA-1A Antibody, Integrin gene promoter Antibody, Alpha polypeptide) Antibody, CD11a antigen Antibody, LFA-1 Antibody, LFA1A Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGAL antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-25
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LFA-1 / CD11A antibody was raised against Human
Antigen Type:	Cells
Immunogen:	LFA-1 / CD11A antibody was raised against leukocytes from a patient suffering from a LGL-type leukaemia.
Specificity:	Reacts with CD11a (alpha-subunit of human LFA-1), a 170-180 kD type I transmembrane glycoprotein expressed on B and T lymphocytes, monocytes, macrophages, neutrophils, basophils and eosinophils.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014 © 2014 LifeSpan BioSciences	