

ITGAM / CD11b Mouse anti-Human Monoclonal (PE) (MEM-174) Antibody - LS-C46438 - LSBio	
CatalogID:	LS-C46438
Target:	integrin, alpha M (complement component 3 receptor 3 subunit) (ITGAM)
Synonyms:	ITGAM Antibody, Antigen CD11b (p170) Antibody, CD11B Antibody, CD11b antigen Antibody, CR3A Antibody, Integrin alpha-M Antibody, MAC-1 Antibody, MO1A Antibody, Neutrophil adherence receptor Antibody, MAC1A Antibody, SLEB6 Antibody, CR-3 alpha chain Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGAM antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MEM-174
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ITGAM / CD11b antibody was raised against Human
Antigen Type:	Cells
Immunogen:	ITGAM / CD11b antibody was raised against human granulocytes.
Specificity:	Recognizes CD11b antigen (Mac-1), a 165 kD type 1 transmembrane protein mainly expressed on monocytes, granulocytes and NK-cells. The CD11b mediates neutrophil and monocyte interactions with stimulated endothelium.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	