

CD14 Mouse anti-Human Monoclonal (aa57-64) (PE) (MEM-18) Antibody - LS-C46439 - LSBio

CatalogID:	LS-C46439
Target:	CD14 molecule
Synonyms:	CD14 Antibody, CD14 molecule Antibody, CD14 antigen Antibody
Host	CD14 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-18
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD14 antibody was raised against Human
Antigen Type:	Protein extract
Immunogen:	CD14 antibody was raised against a crude mixture of human urinary proteins precipitated by ammonium sulphate from the urine of a patient suffering from proteinuria.
Specificity:	Reacts with CD14, a 53-55 kD GPI (glycosylphosphatidylinositol)-linked membrane glycoprotein expressed on monocytes, macrophages and weakly on granulocytes; also expressed by most tissue macrophages. In human, the epitope recognized by MEM-18 is located between amino acids 57-64.
Epitope:	aa57-64
Reactivity:	Human, Primate
Purification:	Ascites
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences