

CD43 Mouse anti-Human Monoclonal (PE) (MEM-59) Antibody - LS-C46445 - LSBio

CatalogID:	LS-C46445
Target:	sialophorin (SPN)
Synonyms:	SPN Antibody, GPL115 Antibody, Leukosialin Antibody, Leukocyte sialoglycoprotein Antibody, Galactoglycoprotein Antibody, GALGP Antibody, Sialophorin Antibody, LSN Antibody, CD43 Antibody, CD43 antigen Antibody
Host	SPN antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-59
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD43 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD43 antibody was raised against human T lymphocytes.
Specificity:	Recognizes neuraminidase-sensitive epitope on CD43 (Leukosialin), a 95-135 kD type I transmembrane glycoprotein (mucin-type) which is involved in lymphocyte activation. CD43 is expressed by platelets and at high levels on the surface of all leukocytes; it is negative on resting B lymphocytes and erythrocytes.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences