

**CD45 Mouse anti-Human Monoclonal (PE) (MEM-28) Antibody - LS-C46447 - LSBio**

<b>CatalogID:</b>	LS-C46447
<b>Target:</b>	protein tyrosine phosphatase, receptor type, C (PTPRC)
<b>Synonyms:</b>	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
<b>Family / Subfamily:</b>	Protein Phosphatase / R1/R6
<b>Host</b>	PTPRC antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	MEM-28
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD45 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45 antibody was raised against human thymocytes and T lymphocytes.
<b>Specificity:</b>	Reacts with all alternative forms of human CD45 antigen (Leukocyte Common Antigen), a 180-220 kD single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except erythrocytes and platelets.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Ascites
<b>Presentation:</b>	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences