

CD45RB Mouse anti-Human Monoclonal (PE) (MEM-55) Antibody - LS-C46449 - LSBio

CatalogID:	LS-C46449
Target:	CD45RB
Host	CD45RB antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-55
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45RB antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RB antibody was raised against human thymocytes and T lymphocytes.
Specificity:	Recognizes a sialidase-sensitive epitope of CD45RB, a 180-240 kD single chain type I membrane glycoprotein, variant of CD45 (CD45RB isoform). CD45RB is expressed on a subset of T lymphocytes, B lymphocytes, monocytes, macrophages, granulocytes and dendritic cells.
Reactivity:	Human, Primate
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences