

CD53 Mouse anti-Human Monoclonal (PE) (MEM-53) Antibody - LS-C46450 - LSBio

CatalogID:	LS-C46450
Target:	CD53 molecule
Synonyms:	CD53 Antibody, CD53 antigen Antibody, CD53 molecule Antibody, CD53 tetraspan antigen Antibody, Cell surface antigen Antibody, Leukocyte surface antigen CD53 Antibody, Transmembrane glycoprotein Antibody, Tetraspanin-25 Antibody, Tspan-25 Antibody, CD53 glycoprotein Antibody, Cell surface glycoprotein CD53 Antibody, MOX44 Antibody, TSPAN25 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD53 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-53
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD53 antibody was raised against Human
Immunogen:	CD53 antibody was raised against leukocytes of patient suffering from a LGL-type leukemia.
Specificity:	Reacts with CD53, a 32-40 kD tetraspanin family glycoprotein exclusively expressed on leukocytes; it is not present on platelets, red blood cells and non-hematopoietic cells. Reacts also with deglycosylated molecule (molecular weight of the antigen is reduced by 15 kD using endoglycosidase F).
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences