

**CD71 / Transferrin Receptor Mouse anti-Human Monoclonal (PE) (MEM-75) Antibody - LS-C46455
- LSBio**

CatalogID:	LS-C46455
Target:	transferrin receptor (TFRC)
Synonyms:	TFRC Antibody, CD71 Antibody, TFR1 Antibody, TRFR Antibody, TFR Antibody, TR Antibody, CD71 antigen Antibody, p90 Antibody, T9 Antibody, Transferrin receptor Antibody, Transferrin receptor protein 1 Antibody
Family / Subfamily:	Protease / Metallopeptidase M28B
Host	TFRC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-75
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD71 / Transferrin Receptor antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD71 / Transferrin Receptor antibody was raised against nALM-6 human pre-B cell line.
Specificity:	Reacts with CD71 antigen (transferrin receptor), a 95 kD type II homodimeric transmembrane glycoprotein expressed on activated B and lymphocytes, macrophages and erythroid precursors; it is lost on resting blood leukocytes. The antibody MEM-75 does not block binding of transferrin to the receptor.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences