

Human Lambda Light Chain Mouse anti-Human Monoclonal (PE) (4C2) Antibody - LS-C46461 - LSBio

CatalogID:	LS-C46461
Target:	Human Ig Lambda Light Chain
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	4C2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Human
Specificity:	Reacts with lambda light chains (22.5 kD) of human immunoglobulin.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences