

CD22 Mouse anti-Human Monoclonal (PE) (IS7) Antibody - LS-C46462 - LSBio

CatalogID:	LS-C46462
Target:	CD22 molecule
Synonyms:	CD22 Antibody, CD22 molecule Antibody, CD22 antigen Antibody, SIGLEC2 Antibody, SIGLEC-2 Antibody, T-cell surface antigen Leu-14 Antibody, B-cell receptor CD22 Antibody, BL-CAM Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD22 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	IS7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD22 antibody was raised against Human
Immunogen:	CD22 antibody was raised against reh human cell line.
Specificity:	Reacts with CD22 (BL-CAM), a 130 kD type I transmembrane glycoprotein (immunoglobulin superfamily) expressed in the cytoplasm of pro-B and pre-B lymphocytes, and on the surface of mature and activated B lymphocytes; it is lost on plasma cells, peripheral blood T lymphocytes, granulocytes and monocytes.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences