

**CD19 Mouse anti-Human Monoclonal (PE) (LT19) Antibody - LS-C46483 - LSBio**

<b>CatalogID:</b>	LS-C46483
<b>Target:</b>	CD19 molecule
<b>Synonyms:</b>	CD19 Antibody, B-lymphocyte antigen CD19 Antibody, CD19 molecule Antibody, Differentiation antigen CD19 Antibody, T-cell surface antigen Leu-12 Antibody, B4 Antibody, CD19 antigen Antibody, CVID3 Antibody
<b>Host</b>	CD19 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	LT19
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD19 antibody was raised against Human
<b>Immunogen:</b>	CD19 antibody was raised against daudi human Burkitt's lymphoma cell line.
<b>Specificity:</b>	Reacts with CD19 (B4), a 95 kD type I transmembrane glycoprotein (immunoglobulin superfamily) expressed on B lymphocytes and follicular dendritic cells; it is lost on plasma cells.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Ascites
<b>Presentation:</b>	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences