

CD27 Mouse anti-Human Monoclonal (PE) (LT27) Antibody - LS-C46485 - LSBio

CatalogID:	LS-C46485
Target:	CD27 molecule
Synonyms:	CD27 Antibody, CD27 antigen Antibody, CD27 molecule Antibody, S152 Antibody, T14 Antibody, T-cell activation antigen CD27 Antibody, CD27L receptor Antibody, T cell activation antigen S152 Antibody, TNFRSF7 Antibody, Tp55 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD27 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	LT27
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD27 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD27 antibody was raised against human peripheral blood lymphocytes.
Specificity:	Reacts with CD27 (T14), a 50-55 kD type I transmembrane glycoprotein (member of the TNF-receptor superfamily) expressed on medullary thymocytes, peripheral T lymphocytes, some B lymphocytes and NK cells.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences