

ITGA2B / CD41 Mouse anti-Human Monoclonal (PE) (MEM-06) Antibody - LS-C46486 - LSBio	
CatalogID:	LS-C46486
Target:	integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41) (ITGA2B)
Synonyms:	ITGA2B Antibody, BDPLT2 Antibody, CD41B Antibody, GT Antibody, gp2B Antibody, Integrin alpha-IIb Antibody, ITGAB Antibody, GPIIb Antibody, HPA3 Antibody, Platelet-specific antigen BAK Antibody, CD41 Antibody, CD41 antigen Antibody, GPalpha IIb Antibody, GTA Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGA2B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-06
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ITGA2B / CD41 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	ITGA2B / CD41 antibody was raised against leukocytes of patient suffering from LGL-type leukaemia.
Specificity:	Reacts with CD41 (GPIIb), a transmembrane glycoprotein (integrin family) composed of two chains GPIIbalpha (heavy chain; 120 kD) and GPIIbbeta (light chain; 23 kD). CD41 is mainly expressed on platelets and megakaryocytes.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	