

MCP / CD46 Mouse anti-Human Monoclonal (PE) (MEM-258) Antibody - LS-C46495 - LSBio	
CatalogID:	LS-C46495
Target:	CD46 molecule, complement regulatory protein
Synonyms:	CD46 Antibody, AHUS2 Antibody, CD46 antigen Antibody, MCP Antibody, Measles virus receptor Antibody, TLX Antibody, MIC10 Antibody, TRA2.10 Antibody, Membrane cofactor protein Antibody
Host	CD46 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MEM-258
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MCP / CD46 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	MCP / CD46 antibody was raised against hPB-ALL human T cell line.
Specificity:	Recognizes an epitope on SCR4 (the membrane-proximal SCR) domain of CD46 (Membrane cofactor protein). CD46 is 56-66 kD dimeric transmembrane protein expressed on T and B lymphocytes, platelets, monocytes, granulocytes, endothelial cells, epithelial cells and fibroblast; it is negative on erythrocytes.
Reactivity:	Human, Bovine
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	