

CD38 Mouse anti-Human Monoclonal (PE) (HIT2) Antibody - LS-C46504 - LSBio

CatalogID:	LS-C46504
Target:	CD38 molecule
Synonyms:	CD38 Antibody, ADP-ribosyl cyclase 1 Antibody, CD38 antigen (p45) Antibody, CD38 antigen Antibody, Cyclic ADP-ribose hydrolase 1 Antibody, T10 Antibody, CADPr hydrolase 1 Antibody, CD38 molecule Antibody, NAD(+) nucleosidase Antibody
Host	CD38 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	HIT2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD38 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD38 antibody was raised against human thymocytes in fetus.
Specificity:	Reacts with CD38 (T10), a 45 kD type II transmembrane glycoprotein strongly expressed mainly on plasma cells and activated T and B lymphocytes; it is an antigenic marker of lymphoid cells.
Reactivity:	Human
Purification:	Ascites
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences