

ANPEP / CD13 Mouse anti-Human Monoclonal (PE) (WM15) Antibody - LS-C46509 - LSBio

CatalogID:	LS-C46509
Target:	alanyl (membrane) aminopeptidase (ANPEP)
Synonyms:	ANPEP Antibody, Alanyl aminopeptidase Antibody, Aminopeptidase M Antibody, APN Antibody, Aminopeptidase N Antibody, AP-N Antibody, AP-M Antibody, CD13 Antibody, CD13 antigen Antibody, HAPN Antibody, gp150 Antibody, p150 Antibody, PEPN Antibody, Microsomal aminopeptidase Antibody, LAP1 Antibody
Family / Subfamily:	Exopeptidase / Metallopeptidase M1
Host	ANPEP antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	WM15
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ANPEP / CD13 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	ANPEP / CD13 antibody was raised against human AML cells.
Specificity:	Recognizes the human CD13 cell surface glycoprotein, a 150kD molecule expressed on granulocytes, endothelial cells, epithelial cells and myeloid progenitors
Reactivity:	Human, Primate
Purification:	Protein G purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences