

CD16 Mouse anti-Human Monoclonal (PE) (LNK16) Antibody - LS-C46510 - LSBio

CatalogID:	LS-C46510
Target:	CD16
Host	CD16 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	LNK16
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD16 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD16 antibody was raised against normal human peripheral blood granulocytes.
Specificity:	Reacts with CD16, a low affinity receptor for aggregated IgG (FcγRIII antigen). CD16 exists in two different isoforms: CD16a (FcγRIIIA; 50-65 kD; expressed on NK-cells, monocytes and macrophages) and CD16b (FcγRIIIB; 48 kD; mainly expressed on neutrophils).
Reactivity:	Human, Primate
Purification:	Ion exchange chromatography
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences