

CD40 Mouse anti-Human Monoclonal (N-Terminus) (PE) (HI40a) Antibody - LS-C46515 - LSBio	
CatalogID:	LS-C46515
Target:	CD40 molecule, TNF receptor superfamily member 5
Synonyms:	CD40 Antibody, B-cell surface antigen CD40 Antibody, B cell-associated molecule Antibody, B cell surface antigen CD40 Antibody, CD40 antigen Antibody, CD40 type II isoform Antibody, p50 Antibody, Bp50 Antibody, CD40L receptor Antibody, CDW40 Antibody, TNFRSF5 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD40 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	HI40a
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD40 antibody was raised against Human
Specificity:	Recognizes CD40 (BP50), a 48 KD type I single chain transmembrane glycoprotein expressed on normal and neoplastic B cells, but not on terminally differentiated plasma cells. CD40 antigen is also present on Hodgkin's and Reed-Sternberg cells, follicular dendritic cells, some macrophages, basal epithelial cells and endothelial cells
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/1/2014 © 2014 LifeSpan BioSciences	