

CD30 Mouse anti-Human Monoclonal (PE) (MEM-268) Antibody - LS-C46526 - LSBio

CatalogID:	LS-C46526
Target:	tumor necrosis factor receptor superfamily, member 8 (TNFRSF8)
Synonyms:	TNFRSF8 Antibody, CD30L receptor Antibody, CD30 Antibody, Cytokine receptor CD30 Antibody, D1S166E Antibody, Ki-1 Antibody, Ki-1 antigen Antibody, CD30 antigen Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF8 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	MEM-268
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD30 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD30 antibody was raised against expression vector containing CD30 cDNA (booster suspension of THP-1 cell line).
Specificity:	Recognizes extracellular part of CD30 (Ki-1 antigen), a 105 kD single chain glycoprotein expressed on Hodgkin's and Reed-Sternberg cells; it is also found in Burkitt's lymphomas, virus-infected T and B lymphocytes, and on normal B and T lymphocytes after activation (T lymphocytes that produce Th2-type cytokines and on CD4+/CD8+ T lymphocytes that co-express CD45RO and the IL4 receptor).
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 15 mM sodium azide, 0.2% high-grade protease free BSA as a stabilizing agent.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences