

CD19 Mouse anti-Human Monoclonal (PE) Antibody - LS-C4713 - LSBio	
CatalogID:	LS-C4713
Target:	CD19 molecule
Synonyms:	CD19 Antibody, B-lymphocyte antigen CD19 Antibody, CD19 molecule Antibody, Differentiation antigen CD19 Antibody, T-cell surface antigen Leu-12 Antibody, B4 Antibody, CD19 antigen Antibody, CVID3 Antibody
Host	CD19 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD19 antibody was raised against Human
Immunogen:	CD19 antibody was raised against protein.
Specificity:	Recognizes human CD19.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow cytometry, Immunohistochemistry and Immunoprecipitation. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	