

CR2 / CD21 Mouse anti-Human Monoclonal (PE) Antibody - LS-C4745 - LSBio

CatalogID:	LS-C4745
Target:	complement component (3d/Epstein Barr virus) receptor 2 (CR2)
Synonyms:	CR2 Antibody, C3DR Antibody, CD21 antigen Antibody, CVID7 Antibody, Complement receptors type 2 Antibody, CD21 Antibody, Complement C3d receptor Antibody, EBV receptor Antibody, SLEB9 Antibody, Complement receptor type 2 Antibody, CR Antibody, Epstein-Barr virus receptor Antibody
Host	CR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CR2 / CD21 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CR2 / CD21 antibody was raised against human IM9 cell line.
Specificity:	Recognizes the human CD21 cell surface antigen, a 145kD glycoprotein expressed by mature B cells and follicular dendritic cells. CD21 acts as a receptor for complement components C3d, C3dg and iC3b and also for Epstein Barr Virus.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	Lyophilized from PBS, pH 7.2, 0.1% sodium azide, 1% BSA.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Flow Cytometry: Neat-1:10.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences