

CD22 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C4765 - LSBio	
CatalogID:	LS-C4765
Target:	CD22 molecule
Synonyms:	CD22 Antibody, CD22 molecule Antibody, CD22 antigen Antibody, SIGLEC2 Antibody, SIGLEC-2 Antibody, T-cell surface antigen Leu-14 Antibody, B-cell receptor CD22 Antibody, BL-CAM Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD22 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD22 antibody was raised against Mouse
Specificity:	Recognizes CD22 (Mr 150kD), a B lymphocyte adhesion molecule.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Blocking of B cell homotypic adhesion. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	