

His Tag Mouse Monoclonal (FITC) Antibody - LS-C48103 - LSBio

CatalogID:	LS-C48103
Target:	His Tag
Host	His Tag antibody was produced in Mouse
Clonality:	Monoclonal
Conjugations:	Fluorescein (FITC)
Purification:	Affinity purified
Presentation:	10 mM phosphate, 150 mM NaCl, pH 7.2, with BSA and 2 mM sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences