

PECAM-1 / CD31 Rat anti-Mouse Monoclonal (PE) (390) Antibody - LS-C4930 - LSBio

CatalogID:	LS-C4930
Target:	platelet/endothelial cell adhesion molecule 1 (PECAM1)
Synonyms:	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
Host	PECAM1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	390
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PECAM-1 / CD31 antibody was raised against Mouse
Immunogen:	PECAM-1 / CD31 antibody was raised against mouse leukocyte cell line 32D.
Specificity:	Recognizes mouse CD31/PECAM-1, Mr 130-140kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Store at 4° C. Do not freeze. Product should be stored undiluted, is photosensitive and should be protected from light.
Usage Summary:	Suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry (frozen sections) and Blocking of cell adhesion. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences