

CD34 Mouse anti-Human Monoclonal Antibody - LS-C4959 - LSBio	
CatalogID:	LS-C4959
Target:	CD34 molecule
Synonyms:	CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody
Host	CD34 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	CD34 antibody was raised against Human
Specificity:	Recognizes human CD34. Cross-reactivity with other species has not been tested.
Reactivity:	Human
Presentation:	PBS, pH 7.2, 0.02% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Suitable for use in Flow Cytometry, Western Blot, Immunoprecipitation and Immunocytochemistry.
Uses:	ICC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	1 ea
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	