

CR1 / CD35 Mouse anti-Human Monoclonal (PE) Antibody - LS-C4990 - LSBio

CatalogID:	LS-C4990
Target:	complement component (3b/4b) receptor 1 (Knops blood group) (CR1)
Synonyms:	CR1 Antibody, C3b/C4b receptor Antibody, C3-binding protein Antibody, CD35 Antibody, C3BR Antibody, Knops blood group antigen Antibody, C4BR Antibody, CD35 antigen Antibody, Complement receptor type 1 Antibody, KN Antibody
Family / Subfamily:	Protease / Serine S1
Host	CR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CR1 / CD35 antibody was raised against Human
Specificity:	Recognizes human CD35.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	+4°C, avoid freezing
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections and paraffin section). Western blotting. Immunoprecipitation. Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Paraffin, IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences