

CD38 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C5062 - LSBio

CatalogID:	LS-C5062
Target:	CD38 molecule
Synonyms:	CD38 Antibody, ADP-ribosyl cyclase 1 Antibody, CD38 antigen (p45) Antibody, CD38 antigen Antibody, Cyclic ADP-ribose hydrolase 1 Antibody, T10 Antibody, CADPr hydrolase 1 Antibody, CD38 molecule Antibody, NAD(+) nucleosidase Antibody
Host	CD38 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD38 antibody was raised against Mouse
Specificity:	Recognizes CD38, Mr 42-46kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry and Induction of B lymphocyte proliferation. Flow Cytometry: 0.2 ug/1E6 cells. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences