

CD38 Rat anti-Mouse Monoclonal (PE, Cy5) Antibody - LS-C5076 - LSBio

CatalogID:	LS-C5076
Target:	CD38 molecule
Synonyms:	CD38 Antibody, ADP-ribosyl cyclase 1 Antibody, CD38 antigen (p45) Antibody, CD38 antigen Antibody, Cyclic ADP-ribose hydrolase 1 Antibody, T10 Antibody, CADPr hydrolase 1 Antibody, CD38 molecule Antibody, NAD(+) nucleosidase Antibody
Host	CD38 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD38 antibody was raised against Mouse
Specificity:	Recognizes CD38 (Mr 42kD), an immunoregulatory ectoenzyme.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. Immunohistochemistry (frozen sections). Induction of B lymphocyte proliferation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences