

HCV NS4 Mouse anti-Hepatitis C Virus Monoclonal (aa1658-1863) (FITC) Antibody - LS-C50904 - LSBio	
<b>CatalogID:</b>	LS-C50904
<b>Target:</b>	Hepatitis C Virus NS4 protein (NS4)
<b>Synonyms:</b>	HCV NS4 Antibody, NS4 Antibody
<b>Host</b>	antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	HCV NS4 antibody was raised against Hepatitis C Virus
<b>Immunogen:</b>	HCV NS4 antibody was raised against hCV NS4 recombinant antigen NS4a+b aa 1658-1863, 19kD.
<b>Specificity:</b>	Specific for HCV NS4 antigen.
<b>Epitope:</b>	aa1658-1863
<b>Reactivity:</b>	Hepatitis C Virus
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.1% sodium azide.
<b>Recommended Storage:</b>	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
<b>Usage Summary:</b>	Western Blot, ELISA, Immunofluorescence and IHC.
<b>Uses:</b>	IHC, Immunofluorescence, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 10/2/2014 © 2014 LifeSpan BioSciences	