

CD40 Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C5130 - LSBio	
CatalogID:	LS-C5130
Target:	CD40 molecule, TNF receptor superfamily member 5
Synonyms:	CD40 Antibody, B-cell surface antigen CD40 Antibody, B cell-associated molecule Antibody, B cell surface antigen CD40 Antibody, CD40 antigen Antibody, CD40 type II isoform Antibody, p50 Antibody, Bp50 Antibody, CD40L receptor Antibody, CDW40 Antibody, TNFRSF5 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD40 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD40 antibody was raised against Mouse
Specificity:	Recognizes mouse CD40, Mr 45-50kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. In vitro induction of proliferation of mature B lymphocytes. Flow Cytometry: 2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	