

CD40 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C5132 - LSBio

CatalogID:	LS-C5132
Target:	CD40 molecule, TNF receptor superfamily member 5
Synonyms:	CD40 Antibody, B-cell surface antigen CD40 Antibody, B cell-associated molecule Antibody, B cell surface antigen CD40 Antibody, CD40 antigen Antibody, CD40 type II isoform Antibody, p50 Antibody, Bp50 Antibody, CD40L receptor Antibody, CDW40 Antibody, TNFRSF5 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD40 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD40 antibody was raised against Mouse
Specificity:	Recognizes mouse CD40, Mr 45-50kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. In vitro induction of proliferation of mature B lymphocytes. Flow Cytometry: 0.5 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences