

CD45RA Mouse anti-Human Monoclonal (LT45RA) Antibody - LS-C51905 - LSBio

CatalogID:	LS-C51905
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	LT45RA
Immunogen Species:	CD45RA antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RA antibody was raised against normal human blood lymphocytes.
Specificity:	CD45RA antigen (MW 220 kD). Working dilution 1:1000.
Reactivity:	Human
Presentation:	PBS, pH 7.4, 0.1 % sodium azide.
Recommended Storage:	+4°C, avoid freezing
Uses:	Undetermined
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences