

CD55 Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C5358 - LSBio	
CatalogID:	LS-C5358
Target:	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
Synonyms:	CD55 Antibody, CD55 antigen Antibody, Decay-accelerating factor Antibody, TC Antibody, CR Antibody, CROM Antibody, DAF Antibody
Host	CD55 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD55 antibody was raised against Human
Specificity:	Recognizes human CD55.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/2/2014 © 2014 LifeSpan BioSciences	