

NCAM / CD56 Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C5398 - LSBio

CatalogID:	LS-C5398
Target:	neural cell adhesion molecule 1 (NCAM1)
Synonyms:	NCAM1 Antibody, CD56 Antibody, CD56 antigen Antibody, MSK39 Antibody, NCAM-1 Antibody, N-CAM-1 Antibody, NCAM Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	NCAM1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	NCAM / CD56 antibody was raised against Human
Specificity:	Recognizes human CD56.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Flow Cytometry: 5 ul/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences