

CD59 Mouse anti-Human Monoclonal (PE) Antibody - LS-C5431 - LSBio	
CatalogID:	LS-C5431
Target:	CD59 molecule, complement regulatory protein
Synonyms:	CD59 Antibody, 16.3A5 Antibody, 1F5 Antibody, 1F5 antigen Antibody, CD59 glycoprotein Antibody, EL32 Antibody, EJ16 Antibody, G344 Antibody, HRF-20 Antibody, MAC-IP Antibody, MACIF Antibody, Ly-6-like protein Antibody, MEM43 Antibody, MIRL Antibody, MEM43 antigen Antibody, MIN2 Antibody, MIN3 Antibody, MSK21 Antibody, p18-20 Antibody, Protectin Antibody, MAC-inhibitory protein Antibody, MIC11 Antibody, T cell-activating protein Antibody, CD59 antigen Antibody, EJ30 Antibody, HRF20 Antibody, Human leukocyte antigen MIC11 Antibody, Lymphocytic antigen CD59/MEM43 Antibody, MIN1 Antibody
Host	CD59 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD59 antibody was raised against Human
Specificity:	Reacts with a PI-linked glycoprotein, Mw 18-20kD found on all types of leucocytes including platelets. The antigen is expressed on many hemopoietic and non-hemopoietic cells.
Reactivity:	Human
Purification:	Protein A purified
Reconstitution:	deionized water
Presentation:	Lyophilized.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Flow Cytometry: Neat.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	