

CD72 Mouse anti-Mouse Monoclonal (PE) Antibody - LS-C5520 - LSBio

CatalogID:	LS-C5520
Target:	CD72 molecule
Synonyms:	CD72 Antibody, CD72 molecule Antibody, CD72 antigen Antibody, CD72b Antibody, Lyb-2 Antibody, LYB2 Antibody
Host	CD72 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD72 antibody was raised against Mouse
Immunogen:	CD72 antibody was raised against dBA/2 mouse spleen cells.
Specificity:	Recognizes mouse CD72.1/Lyb-2.1, 45kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. In vitro blocking and proliferation assays. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences