

CD79B Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C5547 - LSBio

CatalogID:	LS-C5547
Target:	CD79b molecule, immunoglobulin-associated beta
Synonyms:	CD79B Antibody, AGM6 Antibody, B29 Antibody, Ig-beta Antibody, IGB Antibody, CD79b antigen Antibody
Host	CD79B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD79B antibody was raised against Human
Specificity:	Recognizes human CD79beta.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences