

Streptococcus pneumoniae Rabbit anti-Streptococcus pneumoniae Polyclonal (FITC) Antibody - LS-C55807 - LSBio

CatalogID:	LS-C55807
Target:	Streptococcus pneumoniae
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Streptococcus pneumoniae
Reactivity:	Streptococcus pneumoniae
Purification:	Purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 0.1% BSA.
Uses:	Immunofluorescence (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences