

CD80 Rat anti-Mouse Monoclonal (PE, Cy5) Antibody - LS-C5589 - LSBio

CatalogID:	LS-C5589
Target:	CD80 molecule
Synonyms:	CD80 Antibody, Activation B7-1 antigen Antibody, B7-1 Antibody, B7.1 Antibody, BB1 Antibody, CD28LG Antibody, CD28LG1 Antibody, CD80 antigen Antibody, LAB7 Antibody, CD80 molecule Antibody, CTLA-4 counter-receptor B7.1 Antibody
Host	CD80 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD80 antibody was raised against Mouse
Specificity:	Recognizes mouse CD80 (B7-1), Mr 45-47kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. In vivo and in vitro blocking of ligand binding. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences